
Metenova MixOne

SUA 40L Mixing Carboy Truelev

A closed single-use mixing system intended for mixing of pharmaceutical liquid products. The built-in rigidity of the container ensures safe handling and protection of the product inside making it suitable for applications such as transport and storage. The system is designed to promote optimal mixing performance essential for critical pharma applications. The system is tailor made to meet your specific needs of integration.





Container setup

Container	40L Square PP	
Сар	Sealed system (welded lid)	
Mixing system	Truelev 100	
Inlet port	1x3/8"	
Vent Port	1 x 3/8"	
Sampling port	1x 1/8"	
Outlet port	1x 3/8"	
Extra feature 1	Barb protections	
Extra feature 2	N/A	
Extra feature 3	N/A	

General information

Measurements	Max volume (L)		40
	Package mm (in)	Width	580 (22.8)
		Height	900 (35.4)
		Depth	380 (15)
Mixing System	Mixing Unit		Trulev
	Mixing Volume (L)		2,5 - 40
	Compatible drive unit platform		ZU40L513SQT-05171
Material Compliance	Biological reactivity		All product contact material have passed USP VI testing (USP<88>)
	BSE-TSE status		All product contact materials comply to EU guidance EMA/410/01
	Endotoxin		USP <85> (< 0.25 EU/ml)
	Particulate matter		N.A
	Irradiation		Gamma
Assembly	Inspection		100 % leak tested
	Final assembly		Cleanroom ISO Class 7
	Packaging		Double bags in cardboard box







Tube Setup

Inlet line	Length mm (in):	610 (24)
	Tubing	3/8" ID Quick Seal 24" long All pure C-flex
	Fitting	Aseptiquik G Conn. 3/8" HB (AQG17006)
	Fastener	PP tube clamp, 3/8"-3/4"
Outlet line	Length mm (in):	620 (24,4)
	Tubing	3/8" ID 25" Long C-flex 374
	Fitting	Aseptiquik G Conn. 3/8" HB (AQG17006)
	Fastener	PP tube clamp, 3/8"-3/4"
Vent line	Length mm (in):	610 (24)
	Tubing	3/8" ID 25" Long C-flex 374
	Sensor	N.A
	Filter	Opticap XL 50 Express SPG 0,2μm, HB/HB
Sampling line	Length mm (in):	305 (12)
	Tubing	1/8" ID Quick Seal 12" long All pure C-flex
	Sensor	End plug 1/8"
	Fastener	PP tube clamp, 3/8"-3/4"





